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
Vicon Nexus 1.8.5
File Edit Window Help
VICON NEXUS
View Type: Untitled *
Resources
Connected
Go Offline
System
Subjects
ViconTeamViewer *
Local Vicon System [100Hz]
Vicon Cameras
#1 (Bonita 3)
#2 (Bonita 3)
#3 (Bonita 3)
#4 (Bonita 3)
#5 (Bonita 3)
#6 (Bonita 3)
#7 (Bonita 3)
#8 (Bonita 3)
#9 (Bonita 3)
#10 (Bonita 3)
Vicon Connectivity
#1 (MX Giganet)
Analog Card [1000Hz] (Slot...
Devices
Properties
Show Advanced
3D Perspective
LIVE
Tools
Wand: Active Wand
L-Frame: Active Wand
Video Calibration Setup
Activate
Mask Cameras
Aim Cameras
Calibrate Cameras
Show Advanced
Calibration Type: Full Calibration
Cameras To Calibrate: All Cameras
Refinement frames: 1000
DV Calibration frames: 500
Auto Stop:
Start
Cancel
Set Volume Origin
Static Video Calibration
Manage Camera Calibration
Reset
Load...
Camera Calibration Feedback
0%
Camera
Wand Count
Image Error
EN
02:49 ص
1/12/2021

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

Upper Limb (Right)
RSHO – Right Shoulder
RUPA – Right Upper Arm
RELB – Right Elbow
RFRM – Right Forearm
RWRIST – Right Wrist
RTHUMB – Right Thumb
RPINKY – Right Little Finger
Lower Limb (Right)
RTHI – Right Thigh
RKNE – Right Knee
RSHN – Right Shin
RANK – Right Ankle
RTOE – Right Big Toe
RHEE – Right Heel
RMT5 – Right Little Toe

Upper Limb (Right)
RSHO – Right Shoulder
RUPA – Right Upper Arm
RELB – Right Elbow
RFRM – Right Forearm
RWRIST – Right Wrist
RTHUMB – Right Thumb
RPINKY – Right Little Finger
Head
ARIEL – Top of Head
RFHD – Right Forehead
LFHD – Left Forehead
Torso
LFTShould – Left Front Shoulder
RFTShould – Right Front Shoulder
STRN – Sternum
Lower Limb (Right)
RTHI – Right Thigh
RKNE – Right Knee
RSHN – Right Shin
RANK – Right Ankle
RTOE – Right Big Toe
RMT5 – Right Little Toe
Waist
LFWT – Left Front Waist
RFWT – Right Front Waist

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







Vicon Nexus 1.8.5
File Edit Window Help
VICON NEXUS
View Type: Untitled *
Connected
Go Offline
System
Subjects
player 1 (PlugInGait)
3D Perspective
LIVE
Subject: player 1 (PlugInGait)
Subject Capture
Start
Cancel
Frames Captured: 0
0:00:00(00)
Subject Calibration
Pipeline: Calibrate Labeling Model
Hide Pipeline
Properties
Show Advanced
0%
Labeling Template Builder
Create Segments:
Add Marker to Segment:
Link Segments:
Manage Sticks:
Add Parameters:
سيتم ظهور اللاعب بهذا الشكل
EN
08:05 ص
2/7/2021


Vicon Nexus 1.8.5
File Edit Window Help
VICON NEXUS
Connected
Go Offline
System
Subjects
player 1 (PlugInGait)
LIVE
Subject Capture
Subject Calibration
Calibrate Labeling Model
Labeling Template Builder
Create Segments:
Add Marker to Segment:
Link Segments:
Manage Sticks:
Add Parameters:


Vicon Nexus 1.8.5
File Edit Window Help
VICON NEXUS
player 1 Cal 01
Runs the "Reconstruct" pipeline
player 1 (PlugInGait)
Subject Capture
Subject Calibration
Calibrate Labeling Model
Labeling Template Builder
يتم الضغط على
run pipeline
وذلك لظهور العينة

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


Vicon Nexus 1.8.5
File Edit Window Help
VICON NEXUS
player 1 Cal 01
player 1 (PlugInGait)
Available Operations:
FileIO
Calibration
Processing
Workstation Operations
Current Pipeline:
Static Plug-in Gait
VPI Compatibility Run static gait model (Run static
Processing Static Subject Calibration (Static Subje
Export data to ASCII file (Advanced)
Export data to ASCII file (Advanced)

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